

Using variable ionisation energy time-of-flight mass spectrometry with comprehensive GC×GC to identify isomeric species

Alam, Mohammed; Stark, Christopher; Harrison, Roy

License:

None: All rights reserved

Document Version

Peer reviewed version

Citation for published version (Harvard):

Alam, M, Stark, C & Harrison, R 2016, 'Using variable ionisation energy time-of-flight mass spectrometry with comprehensive GC×GC to identify isomeric species', *Analytical Chemistry*.

[Link to publication on Research at Birmingham portal](#)

Publisher Rights Statement:

Checked for eligibility: 08/04/2016

General rights

Unless a licence is specified above, all rights (including copyright and moral rights) in this document are retained by the authors and/or the copyright holders. The express permission of the copyright holder must be obtained for any use of this material other than for purposes permitted by law.

- Users may freely distribute the URL that is used to identify this publication.
- Users may download and/or print one copy of the publication from the University of Birmingham research portal for the purpose of private study or non-commercial research.
- User may use extracts from the document in line with the concept of 'fair dealing' under the Copyright, Designs and Patents Act 1988 (?)
- Users may not further distribute the material nor use it for the purposes of commercial gain.

Where a licence is displayed above, please note the terms and conditions of the licence govern your use of this document.

When citing, please reference the published version.

Take down policy

While the University of Birmingham exercises care and attention in making items available there are rare occasions when an item has been uploaded in error or has been deemed to be commercially or otherwise sensitive.

If you believe that this is the case for this document, please contact UBIRA@lists.bham.ac.uk providing details and we will remove access to the work immediately and investigate.

Using variable ionisation energy time-of-flight mass spectrometry with comprehensive GC×GC to identify isomeric species

Mohammed Salim Alam, Christopher Stark, and Roy Michael Harrison

Anal. Chem., **Just Accepted Manuscript** • DOI: 10.1021/acs.analchem.5b03122 • Publication Date (Web): 16 Mar 2016

Downloaded from <http://pubs.acs.org> on March 21, 2016

Just Accepted

“Just Accepted” manuscripts have been peer-reviewed and accepted for publication. They are posted online prior to technical editing, formatting for publication and author proofing. The American Chemical Society provides “Just Accepted” as a free service to the research community to expedite the dissemination of scientific material as soon as possible after acceptance. “Just Accepted” manuscripts appear in full in PDF format accompanied by an HTML abstract. “Just Accepted” manuscripts have been fully peer reviewed, but should not be considered the official version of record. They are accessible to all readers and citable by the Digital Object Identifier (DOI®). “Just Accepted” is an optional service offered to authors. Therefore, the “Just Accepted” Web site may not include all articles that will be published in the journal. After a manuscript is technically edited and formatted, it will be removed from the “Just Accepted” Web site and published as an ASAP article. Note that technical editing may introduce minor changes to the manuscript text and/or graphics which could affect content, and all legal disclaimers and ethical guidelines that apply to the journal pertain. ACS cannot be held responsible for errors or consequences arising from the use of information contained in these “Just Accepted” manuscripts.



1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53

Using variable ionisation energy time-of-flight mass spectrometry with comprehensive GC×GC to identify isomeric species

54
55
56
57
58
59
60

Mohammed S. Alam, Christopher Stark, and
Roy M. Harrison^{*†}

Division of Environmental Health and Risk Management
School of Geography, Earth and Environmental Sciences
University of Birmingham
Edgbaston, Birmingham B15 2TT
United Kingdom

* To whom correspondence should be addressed.

Tele: +44 121 414 3494; Fax: +44 121 414 3708; Email: r.m.harrison@bham.ac.uk

†Also at: Department of Environmental Sciences / Center of Excellence in Environmental Studies, King Abdulaziz University, PO Box 80203, Jeddah, 21589, Saudi Arabia

26 **ABSTRACT**

27 Although GC×GC-ToF-MS allows the separation of thousands of peaks, many of these peaks are
28 not positively identified owing to the lack of mass spectral library data and/or standard materials,
29 leading to a substantial amount of information being inaccessible. The fragmentation patterns of
30 molecules in mass spectrometers using electron impact ionisation at 70 eV can be useful for
31 molecule identification, provided a match is available in a published EI MS library, but are
32 indistinguishable for many isomeric organic compounds (for example, aliphatic and branched
33 alkanes). Lower ionisation energies have been exploited leading to organic compounds being
34 ionised with lower excess internal energy and less fragmentation, retaining the molecular ion and
35 maximising its relative signal. This has enabled the identification of a large number of isomeric
36 organic compounds, both aliphatic and aromatic, between C₁₂-C₃₆, in the previously unresolved
37 complex mixture (UCM) of two motor oil samples.

38
39 This technique also demonstrates problems associated with separation of co-eluting isomers,
40 particularly for the *n*-alkanes, which are routinely measured by 1D GC/MS and may be
41 overestimated, due to co-elution. As a consequence retention times in 2 dimensions and mass
42 spectra at variable ionisation energies are shown to give unparalleled power to identify specific
43 isomers.

44 **Keywords:** GC×GC, gas chromatography, mass spectrometry, ionisation, *n*-alkanes, fragmentation

45
46

47 **INTRODUCTION**

48 There is a well-established need for the characterisation of organic matter in a range of media. The
49 study of organic compounds in the atmosphere has been an interest of scientists for many years
50 owing to the potential detrimental health effects.¹⁻⁴ The complexity of the atmosphere as well as the
51 drawbacks associated with current analytical techniques have limited studies in investigating
52 specific compound classes and/or target analytes. Thus, detailed knowledge of the identities and
53 chemistry of organic compounds as a whole remains elusive. Goldstein and Galbally⁵ illustrated the
54 exponential increase of number of unique isomers possible as a function of the number of carbon
55 atoms in an organic molecule, with *ca.* 1×10^6 constitutional isomers possible for a C₂₀ alkane.

56
57 The study of organic compounds within crude oils has also generated much interest as it plays an
58 important role in petroleum exploration. Crude oils have a limited number of compound classes, but
59 an enormous number of individual components within a class. The composition of crude oil is also
60 susceptible to physical, chemical, biological and weathering processes, leading to significant
61 differences between any two crude oil samples.⁶ Sutton et al.⁷ isolated individual hydrocarbon
62 fractions from an unresolved complex mixture (UCM) of hydrocarbons isolated from a biodegraded
63 crude oil and suggested that the entire oil UCM may contain *ca.* 250,000 unidentified compounds.
64 This suggests that a substantial fraction of chemical composition information is inaccessible and
65 therefore unexploited.

66
67 The most common method for the speciation of organic compounds has been gas chromatography–
68 mass spectrometry (GC/MS). This technique has been exploited to separate crude oils;⁷⁻⁹ airborne
69 samples which encompass multiple emission sources containing multiple compound classes and
70 both polar and non-polar components;¹⁰⁻¹³ and sediments which can contain fingerprints of both
71 crude oils and various other source emissions.^{6,14} The peak capacity of conventional GC, however,
72 is insufficient to resolve the full complexity of these mixtures and a pronounced rising baseline or

1
2 73 series of rising baselines in the gas chromatogram are observed, often referred to as the unresolved
3
4 74 complex mixture (UCM).¹⁵ For example, the UCM in conventional GC separation of oil-based
5
6 75 samples typically accounts for *ca.* 95 % of the area of the chromatogram and is usually unresolved
7
8 76 by analysis. It is thought to comprise a major fraction of constitutional linear, branched and cyclic
9
10 77 alkane isomers and a minor fraction of aromatic compounds.^{16,17} Due to the difficulties in
11
12 78 identifying isomers, information regarding their molecular structures has generally been
13
14 79 unavailable. This knowledge, however, is crucial in understanding the environmental fates of these
15
16 80 compounds, as the degree and positioning of branching in these alkanes strongly affect their
17
18 81 oxidation mechanisms in the environment and secondary organic aerosol (SOA) formation yields in
19
20 82 the atmosphere.¹⁸⁻²⁰
21
22
23

24 83

25
26 84 Multidimensional or comprehensive gas chromatography (GC×GC) provides enhanced peak
27
28 85 capacity and separating powers, due to an orthogonal separation. This technique utilises two
29
30 86 capillary columns of different stationary phases connected in series by a modulator, which delivers
31
32 87 discrete bundles of effluent from the primary column to the secondary column.^{21,22} Compounds of
33
34 88 similar chemical structure are grouped into distinct patterns in ordered chromatograms, providing
35
36 89 useful information that aids identification.²³ GC×GC is often coupled with time-of-flight mass
37
38 90 spectrometry (ToF-MS) which provides extra information in the identification of compounds, due
39
40 91 to its fast data acquisition rate, full range mass sensitivity and, in some cases, high resolution. These
41
42 92 advantages make GC×GC-ToF-MS an attractive analytical technique for separating complex
43
44 93 mixtures, measuring a wide range of volatile and semi-volatile organic compounds in air and oils.²⁴⁻
45
46
47

48 94 ²⁷

49 95

50
51
52 96 Although GC×GC-ToF-MS allows the separation of thousands of peaks, many of these peaks are
53
54 97 not positively identified owing to the lack of mass spectral library data and/or standard materials.²⁸
55
56

57 98 This means that a substantial amount of information is still inaccessible. Traditional mass
58
59
60

1
2 99 spectrometers (MS) employ electron impact (EI) ionisation at 70 eV, where a molecule undergoes
3
4 100 ionisation by an electron (most organic molecules have ionisation energies <11 eV). The electron
5
6 101 imparts a large amount of excess energy causing extensive fragmentation. The fragmentation mass
7
8 102 spectra can in most cases be useful for molecule identification, provided a match is available in a
9
10 103 published EI MS library, but for many organic compounds, such as the isomeric alkane structures
11
12 104 found in the UCM, are indistinguishable. This is because the fragment ions found in the mass
13
14 105 spectra of aliphatic compounds, for example alkanes, are very similar *i.e.* m/z 43, 57, 71, 85, 99 *etc*
15
16 106 and the molecular ion is often barely visible or absent. Interpretation of EI-MS spectra of complex
17
18 107 mixtures can therefore be ambiguous. These limitations in EI have led to a continual demand for
19
20 108 more robust lower (or soft) ionisation techniques for MS. Using lower ionisation energies can lead
21
22 109 to organic compounds being ionised with minimal excess internal energy and thus less
23
24 110 fragmentation, hence retaining the identity of the molecule by retaining a much larger fraction of
25
26 111 the molecular ion. This has been the subject of rigorous research by Maccoll and co-workers,²⁹
27
28 112 since their early work on ion enthalpies and their application in mass spectrometry.³⁰ A series of
29
30 113 literature has since been published, in seventeen parts over a decade, reporting low EI (12.1 eV),
31
32 114 low temperature (~350K) mass spectra of various compound classes.
33
34
35
36
37
38

39 115
40 116 Several soft ionisation techniques have been developed such as chemical ionisation (CI),^{31,32} field
41
42 117 ionisation (FI)^{33,34} and photoionisation (PI).^{35,36} CI typically utilises reagent ions produced by EI of
43
44 118 a reagent gas (such as methane, isobutane or ammonia) at high pressure (*ca.* 1×10^{-4} mbar). The
45
46 119 reagent ions are only slightly reactive with the reagent gas itself, but readily react to ionise the
47
48 120 sample via ion-molecule reactions. Although this ionisation technique often gives molecular mass
49
50 121 information with reduced fragmentation compared to EI, results are highly dependent upon reagent
51
52 122 gas type, pressure, reaction time and the nature of the sample. In addition, CI requires a different
53
54 123 ion-source configuration and thus it can be time consuming and resource intensive if a single MS
55
56 124 instrument is being used for both EI and CI. FI has been exploited by a number of studies.³⁷⁻³⁹ FI
57
58
59
60

1
2 125 produces dominant molecular ions with little or no fragmentation. Hejazi et al.³⁸ developed a
3
4 126 method that utilises GC with parallel hard (EI) and soft ionisation (FI) mass spectrometry. This
5
6 127 instrument consisted of two mass spectrometers, with differing ionisation methods, connected to a
7
8 128 single GC. Wang et al.³⁹ investigated the separation of diesel fuel using a two dimensional
9
10 129 separation approach whereby FI-MS was coupled to a GC. FI, however, produces weak ion currents
11
12 130 with significantly more fluctuation on the one second time scale of a typical scan for a scanning
13
14 131 MS, resulting in poor repeatability of the relative abundances of measured ions.³⁷
15
16
17 132

18
19 133 An increasingly popular soft ionisation technique is PI, where principally two different PI
20
21 134 techniques have been investigated, namely resonance-enhanced multi-photon ionisation (REMPI)
22
23 135 and single photon ionisation (SPI).³⁵ REMPI utilises intense ultraviolet (UV) light laser pulses for a
24
25 136 two photon absorption process, whereas SPI uses vacuum ultraviolet (VUV) photons for a single
26
27 137 photon absorption process. REMPI is highly sensitive and selective for specific compounds such as
28
29 138 aromatics and polyaromatic compounds,^{36,40,41} which prevents this soft ionisation technique from
30
31 139 being applicable to crude oil samples that are rich in aliphatic compounds. SPI, however, with VUV
32
33 140 is capable of soft ionisation of all organic compound classes including the aliphatic alkanes.^{35,42}
34
35 141 Isaacman et al.⁴³ used VUV ionisation at 10 – 10.5 eV to elucidate some of the structures within the
36
37 142 UCM of atmospheric particles, by separating components using GC. Chan et al.⁴⁴ also exploited
38
39 143 this technique in order to characterise the UCM and to provide insights into atmospheric processing
40
41 144 of these components. Their results showed that the most abundant isomer for each carbon number
42
43 145 was the linear alkane accounting for >30% of the total alkane mass. Drawbacks of SPI VUV,
44
45 146 however, are the relative experimental complexity of the technique and the reduced stability of the
46
47 147 molecular ions that are formed (radical cations) in comparison to other techniques such as MALDI
48
49 148 and atmospheric pressure CI (even-electron ions). It is therefore, important to develop and
50
51 149 investigate reliable, inexpensive and simple ionisation techniques that can help identify the highly
52
53 150 numerous unknown compounds that are present in samples such as aerosols and petroleum samples.
54
55
56
57
58
59
60

1
2 151 In this study, we use variable ionisation time-of-flight mass spectrometry, ranging from 10 to 70 eV
3
4 152 coupled to GC×GC, to determine the structures of aliphatic and aromatic compounds within the
5
6 153 UCM. By using soft EI and conventional 70 eV EI, both molecular ions and fragment ions can be
7
8 154 obtained. EI provides high ion yield and sensitivity, generating highly reproducible spectra, for
9
10 155 which there is a large database of spectra of known compounds (NIST). Soft ionisation provides the
11
12 156 molecular ion and reduced fragmentation, for which there is no database available. Both mass
13
14 157 spectral data at 14 and 70 eV, as well as retention times from the GC×GC, give insights into
15
16 158 structural positioning and extent of alkyl branching of aliphatic (linear) alkanes and aromatics.
17
18
19
20 159

21 160 **EXPERIMENTAL**

22 161 **Instrumentation**

23
24 162 Sample analyses were conducted by a gas chromatograph (7890B, Agilent Technologies,
25
26 163 Wilmington, DE, USA) equipped with a Zoex ZX2 cryogenic modulator (Houston, TX, USA). The
27
28 164 first dimension was equipped with a SGE DBX5, non-polar capillary column (30 m, 0.25 mm ID,
29
30 165 0.25 μm – 5% phenyl polysilphenylene-siloxane), and the second dimension column was SGE
31
32 166 DBX50 (4.0 m, 0.1 mm ID, 0.1 μm – 50% phenyl polysilphenylene-siloxane). The GC×GC was
33
34 167 interfaced with a BenchTOF-Select, time-of-flight mass spectrometer (Markes International,
35
36 168 Llantrisant, UK). The acquisition speed was 50 Hz with a mass resolution of >1200 FWHM at 70
37
38 169 eV and >800 FWHM at 14eV over 100 – 1000 m/z .⁴⁵ The mass range was 35 to 600 m/z . The GC
39
40 170 sensitivity at 70 eV has been calculated by the injection of 1 pg octafluoronaphthalene (m/z 272)
41
42 171 giving a signal-to-noise ratio (S/N) of >2000:1 rms, while for 14eV, the injection of 1 pg of
43
44 172 benzophenone (m/z 182) gave a S/N of >80:1 rms.⁴⁵ Although ionisation energies can be varied
45
46 173 from anywhere between 10 and 70 eV, explored energies in this work were 14 eV and 70 eV, by
47
48 174 means of repeat injection runs. All data produced were processed using GC Image v2.3 (Zoex
49
50 175 Corporation, Houston, US).
51
52
53
54
55
56
57
58
59
60

1
2 177 **Principle of soft ionisation technique**

3
4 178 For the conventional EI ion source, electrons are accelerated from the surface of a negatively
5
6 179 charged filament towards a positively charged ion chamber, owing to a potential difference of 70 V,
7
8 180 thus efficiently ionising compounds from the GC. Reducing this potential difference (for soft
9
10 181 ionisation) results in a dramatic loss of sensitivity due to the inefficiency of drawing electrons into
11
12 182 the ion chamber and clustering of electrons around the filament. The variable ionisation source uses
13
14 183 ion optics in order to retain the high potential difference, thus accelerating the electrons from the
15
16 184 filament and reducing their energy prior to arriving at the ion chamber.⁴⁶ This allows the ionisation
17
18 185 energy to be tuned between 10 and 70 eV. Unlike CI and FI, no reagent gases, adjustments in
19
20 186 pressure or switching of sources are required.
21
22
23

24 187
25
26 188 **Methodology**

27
28 189 Two motor/lubricating oil samples (synthetic 5W30 and pure base oil analogous to mineral oil, see
29
30 190 Supporting Information for more details) were analysed, utilising GC×GC-ToF-MS to identify
31
32 191 aliphatic isomeric compounds. The protocol applied to the motor oil analytical sample was to inject
33
34 192 1 µL of diluted (1:1000) sample in a split ratio 100:1 at 300 °C. The initial temperature of the
35
36 193 primary oven (175°C) was held for 5 min and then increased at 1°C/min to 325°C. The secondary
37
38 194 oven temperature program was the same as the primary oven but was increased at 1°C/min to
39
40 195 330°C to ensure all species passed through the column. The transfer line temperature was 325°C
41
42 196 and the ion source temperature was 280°C. Helium was used as the carrier gas at a constant flow
43
44 197 rate of 0.8 mL/min. A modulation time of 7 seconds was used, while a total run time of each sample
45
46 198 was 120 minutes.
47
48
49

50 199
51
52 200 **RESULTS & DISCUSSION**

53
54 201 A typical two-dimensional separation is presented in Figure 1(A) and (B) for the base lubricating
55
56 202 oil. The white line in Figure 1(A) illustrates the separation achievable using traditional 1D GC/MS.
57
58
59
60

1
2 203 Each coloured spot represents an individual compound with a corresponding full mass spectrum.
3
4 204 The intense dark blue spots are major peaks, while the lighter blue spots represent small peaks.
5
6 205 Compounds were separated in the first dimension (X-axis) on a volatility basis and polarity in the
7
8 206 second (Y-axis).
9
10 207
11
12 208 Samples were run using two ionisation energies, 70eV and 14eV to gain insight on fragmentation
13
14 209 patterns and molecular ions, respectively. The procedure for peak identification is outlined in detail
15
16 210 elsewhere.²³ Briefly, peaks were identified based primarily on mass spectral data from the NIST
17
18 211 library (for 70eV), where two match factors (MF>750 and RMF>800) were used to assess the
19
20 212 similarity between the library and measured mass spectra. For samples analysed with lower
21
22 213 ionisation energy (14eV), the identification of species was determined by interpretation of mass
23
24 214 spectral data and comparing their retention indices and mass spectra at 70eV, together with the
25
26 215 presence of the species molecular ion (possessing higher relative intensity than in a 70eV mass
27
28 216 spectrum). The samples analysed at both ionisation energies therefore provided greater insight into
29
30 217 a compound's identity, indicating the power of soft EI. Many differences were observed between
31
32 218 the two oils (total ion chromatogram for 5W30 oil not shown in Figure 1A), the most notable being
33
34 219 the presence of more high molecular weight (HMW) (less volatile) compounds in the 5W30
35
36 220 synthetic oil.

221

222 *Monomethylalkanes*

223 Monomethylalkanes are present in many natural and synthetic organic materials, including crude
224 oil, coal and sediments;⁴⁷⁻⁵⁰ and have been proposed as endogenous biomarkers in exhaled breath.⁵¹
225 The separation and identification of monomethylalkanes in a broad range of carbon atoms is
226 difficult owing to the number of isomeric possibilities and close retention times. Recently, Szabó et
227 al⁵² utilised GC×GC-FID and GC/MS/MS with 100m columns to separate 63 C₉-C₁₉
228 monomethylalkanes, reporting up to 8 isomers for the monomethyl-nonadecane (C₁₉) in exhaled

1
2 229 breath. Krkošová et al.⁴⁷ reported the linear retention indices of 196 C₄-C₃₀ monomethylalkanes
3
4 230 using GC/MS equipped with a 100m column and used these retention indices to identify these
5
6 231 species in diesel fuel. All these studies have considered a targeted approach in their analyses and
7
8 232 not considered the presence of other compounds in complex mixtures.
9

10
11 233

12
13 234 The use of variable ionisation energy coupled to mass spectrometry, as well as retention time
14
15 235 information enables a non-targeted approach to be conducted. For example, 2-methyl-eicosane and
16
17 236 3-methyl-eicosane are C₂₁H₄₄ isomers which have different retention indices and are largely
18
19 237 experimentally indistinguishable from other aliphatic isomers using their fragmentation mass
20
21 238 spectra at 70eV (see Figure 2A and 2D). The absence of the molecular ion (*m/z* 296) means that it is
22
23 239 difficult to deduce, without retention time information, which isomeric aliphatic (with
24
25 240 corresponding C number) it belongs to; as all aliphatic compounds with ≥C₂₂ show relatively small
26
27
28 241 peaks at *m/z* 253 and *m/z* 267. At 14eV however, both these compounds are distinguishable,
29
30 242 allowing the identification of the positioning of branching on an aliphatic alkane (see Figure 2C and
31
32 243 2F). The presence of the molecular ion at *m/z* 296 allows the identification of these compounds
33
34 244 without the use of retention indices. The 14eV mass spectrum of 2-methyl-eicosane shows a
35
36 245 relatively large peak at *m/z* 253, signifying a loss of a C₃H₇ fragment (see Figure 2C), while for 3-
37
38 246 methyl-eicosane the presence of a significant peak at *m/z* 267 signifies the loss of C₂H₅ indicating
39
40 247 the most likely point of cleavage (see Figure 2F). This information together with their identical
41
42 248 retention indices in both 14eV chromatography and 70eV chromatography allows the identification
43
44 249 of 7 further C₂₁ monomethylalkane isomers. For example, the 14eV mass spectrum of 4-methyl-
45
46 250 eicosane showed the presence of peaks *m/z* 225 signifying the loss of a C₅H₁₁ fragment and *m/z* 253,
47
48 251 signifying a loss of a C₃H₇ fragment (see Figure 2J). 2-Methyl and 4-methyl-eicosane mass spectra
49
50 252 were distinguishable from one another, due to the presence of the *m/z* 225 fragment, together with
51
52 253 their positioning in the chromatography which was confirmed using Kovats retention indices. The
53
54 254 mass spectra of 5-methyl-eicosane shows relatively large peaks at *m/z* 239 and *m/z* 211 signifying
55
56
57
58
59
60

1
2 255 the loss of mass fragments C_4H_9 and C_6H_{13} , respectively (see Figure 2K). 6-Methyl-eicosane mass
3
4 256 spectra exhibited peaks at m/z 225 and m/z 197, where the m/z 71 and m/z 99 peaks showed similar
5
6 257 relative intensities (see Figure 2L). From these peaks together with the m/z 197, 6-methyl-eicosane
7
8 258 was distinguishable from 4-methyl-eicosane (which has the same m/z 225 fragment). Peaks at m/z
9
10 259 225 and 197 corresponded to C_5H_{11} and C_7H_{15} fragment losses. This analogy was completed for 7,
11
12 260 8, 9 and 10-methyl-eicosanes; however, there is significant co-elution for 9-methyleicosane and 10-
13
14 261 methyleicosane compounds (Figure 3(A)). All methyl-eicosane mass spectra retained the presence
15
16 262 of the molecular ion at m/z 296 (see Figure 2A-2L).
17
18
19
20 263

21 264 As compounds with similar chemical structure are grouped into distinct patterns in ordered
22
23 265 chromatograms, the order of branching with increasing retention time is consistent throughout the
24
25 266 chromatogram with increasing carbon number (up to C_{28}). The faster eluting monomethylalkanes
26
27 267 are those methylated alkanes with methyl groups furthest away from the first carbon position (with
28
29 268 the exception of 2-methylmonoalkanes, determined by assessing each monomethylalkane's mass
30
31 269 spectrum at 14eV and previously published literature.⁴⁷ Thus, for C_{21} isomers the order of branching
32
33 270 with increasing retention time would be 10-, 9-, 8-, 7-, 6-, 5-, 4-, 2- and then 3-methyl-eicosane.
34
35 271 The majority of these isomeric monomethylalkanes are present in the synthetic oil and are not all
36
37 272 present within the base oil. For example, for C_{21} isomers, there is no 5-methyl-eicosane present in
38
39 273 the base oil (see Figure 3A). Signal intensities of the monomethyl C_{24} – C_{28} alkanes in the base oil
40
41 274 are also 20-30% lower than the synthetic oil, which may reflect the differences in the production
42
43 275 and refining of the two oils (the base oil being an API group 1 oil and having ~77% saturates).
44
45 276 These compounds were easily identifiable using 14eV spectra throughout the UCM for the base oil
46
47 277 but not for the 5W30 oil, most probably due to the coelution of the larger number of isomers present
48
49 278 in the latter. The base oil also consists of a number of di- tri- tetra- and penta-methylalkanes
50
51 279 extending beyond each *n*-alkane species as shown in Figure 3(B), identified tentatively by their
52
53 280 mass spectra at both ionisation energies (for $<C_{22}$) and published retention times (where available).
54
55
56
57
58
59
60

1
2 281 The tri- and tetra-methyl alkanes (belonging to the C₂₃ carbon number) were not present in
3
4 282 abundance in the 5W30 lubricating oil (see Figure 3) and elute faster than the C₂₁ *n*-alkane. It must
5
6 283 be noted that the first dimension is separated by volatility and is not a carbon number distribution.
7
8 284 The mass spectra for many substituted methylalkanes with >C₂₂ were indistinguishable at 70 eV.
9
10 285 Figure 4 illustrates the experimentally determined mass spectra for 2,6,10,14-tetramethyl-
11
12 286 nonadecane (C₂₃H₄₈; upper panel) and 2,6,10,14,18-pentamethyl-nonadecane (C₂₄H₅₀; lower panel)
13
14 287 at 70 eV (4A), 14 eV (4B) and their respective NIST reference spectra (4C). Although these
15
16 288 compounds, in both motor oils, elute at their specific published retention times, they are
17
18 289 indistinguishable at 70 eV. The mass spectral peaks are those found in all aliphatic compounds with
19
20 290 several CH₂ (*m/z* 14) cleaved mass fragments. At 14 eV, however, the presence of the molecular ion
21
22 291 and relatively larger intensities of ions that signify likelihood of the cleavage of a C-C bond, allows
23
24 292 the positive identification of these compounds within complex mixtures. For example, for
25
26 293 2,6,10,14-tetramethyl-nonadecane, the loss of fragments C₃H₁₁ + C₆H₁₃, C₁₀H₂₁ + C₁₁H₂₃ and
27
28 294 C₁₅H₃₁ + C₁₆H₃₃ signify bond breaking either side of the methyl group, indicating the positioning of
29
30 295 the methyl carbon bond and give rise to *m/z* 253+239, 183+169, and 113+99, respectively. While
31
32 296 for 2,6,10,14,18-pentamethyl-nonadecane the loss of C₁₁H₂₃ and C₈H₁₇ fragments results in
33
34 297 relatively large peaks at *m/z* 183 and 113, respectively. This information as well as the presence of
35
36 298 the molecular ion enables confident structural assignment of these compounds. It is interesting that
37
38 299 the NIST (70eV) reference spectrum looks significantly softer than the experimental 70eV spectrum
39
40 300 and remarkably similar to the 14eV ionisation mass spectrum. The NIST library at 70 eV also
41
42 301 signifies the presence of the molecular ion (see Figure 4C); however, this is not easily observed for
43
44 302 many aliphatics within complex mixtures, particularly for high molecular weights.⁵³
45
46
47
48
49
50
51
52

53 304 *n*-Alkanes

54
55 305 *n*-Alkanes are ubiquitous in the environment. They are biosynthesized by several organisms such as
56
57 306 terrestrial and aquatic plants, algae and insects and are used as biomarkers in biogeochemistry.⁵⁴ *n*-
58
59
60

1
2 307 Alkanes are a major constituent of crude oil and are routinely and easily measured by GC/MS. In
3
4 308 petroleum polluted sediments the *n*-alkanes generally range between C₁₂ and C₄₀, whereas in marine
5
6 309 organisms the range is normally <C₂₁.⁵⁵
7

8
9 310

10 311 Diagnostic ratios such as, natural *n*-alkane ratio (NAR), carbon preference index (CPI),
11
12 312 terrigenous/aquatic ratio (TAR) and low molecular weight (LMW)/HMW ratios are exploited to
13
14 313 assess possible sources of *n*-alkanes. For example, NAR estimates the proportion of natural and
15
16 314 petroleum *n*-alkanes, where ratios close to zero signify petroleum hydrocarbons and crude oil
17
18 315 sources, while ratios close to one suggest terrestrial or marine plant sources.⁵⁶ CPI are ratios
19
20 316 between odd and even numbered *n*-alkanes, where values in the range *n*-C₂₄ and *n*-C₃₂ close to one
21
22 317 are found for crude oils and other petroleum sources, while a predominance of odd numbered *n*-
23
24 318 alkanes is characteristic of terrestrial plant debris.^{55,57,58} TAR is used as a crude indicator of relative
25
26 319 terrigenous versus aquatic organic matter input, where the ratio between the concentrations of long
27
28 320 chain *n*-alkanes (*n*-C₂₇ + *n*-C₂₉ + *n*-C₃₁) and short chain *n*-alkanes (*n*-C₁₅ + *n*-C₁₇ + *n*-C₁₉) is
29
30 321 evaluated.⁵⁹ All these ratios are crude indicators of the sources of *n*-alkanes where concentrations
31
32 322 are generally measured by conventional GC/MS.
33
34
35
36

37 323

38
39 324 When these hydrocarbons are measured from within a complex mixture using conventional GC/MS,
40
41 325 large distinct peaks are observed above a rising baseline (or UCM). Although the retention times of
42
43 326 *n*-alkanes are very different from one another, their retention times may be similar to other
44
45 327 structural isomers, particularly for larger carbon numbers, owing to the increased probability of
46
47 328 possible isomers. If co-elution occurs, however, then the MS used in conjunction with GC could
48
49 329 potentially distinguish between the two co-eluting peaks, unless they have similar mass spectra at
50
51 330 70eV.
52

53
54
55 331
56
57
58
59
60

1
2 332 Figure 5A and 5B shows the contour plots of two portions of the motor oils, focusing on the *n*-
3
4 333 alkanes between C₁₈ and C₂₂ (see Figure S1 for contour plot for 5W30 for *n*-alkanes C₂₃ to C₃₀).
5
6 334 Figure 5A (synthetic 5W30) shows the presence of an unknown species (compound X) eluting at
7
8 335 the same retention time as the *n*-alkane (in the X-axis), and is separated by polarity in the Y-axis.
9
10 336 The MS of compound X (see Figure 5A) possesses the molecular ion of the C_{*n*} *n*-alkane + CH₂ (*m/z*
11
12 337 14) at 14eV, suggesting that it is a branched alkane belonging to the C_{*n*+1} aliphatic alkane. This is
13
14 338 illustrated by Figure 5A where the black lines signify the overlap between the isomers of two *n*-
15
16 339 alkanes. The mass spectrum of compound X at 70eV, however, cannot be used to distinguish
17
18 340 between the C_{*n*} *n*-alkane and the C_{*n*+1}-alkane isomer, as there is no molecular ion (particularly for
19
20 341 the larger *n*-alkanes, >C₁₉), see Figure 5(A). This unknown species is not distinguishable from the
21
22 342 *n*-alkane (for >C₁₉) when analysing using conventional GC (or by GC×GC with shorter secondary
23
24 343 columns), as the retention time (volatility axis) of compound X is that of the *n*-alkane and therefore
25
26 344 it will co-elute. This co-elution cannot be separated by its MS at 70eV as the primary and secondary
27
28 345 ions that would be used to identify both compound X and the *n*-alkane would be the same (*i.e.* *m/z*
29
30 346 43, 57, 71, 85) see embedded MS in Figure 5A. In Figure 5A, the volume of compound X co-
31
32 347 eluting with the *n*-alkane is approximately 40% of the total volume of the two together.
33
34 348 Furthermore, the separation of these two compounds by polarity (y-axis) is enhanced when
35
36 349 increasing the length of the secondary column (to 4m), as when using a shorter length column, GC
37
38 350 Imagev2.3 software was unable to automatically deconvolve the two peaks. This suggests that in
39
40 351 studies utilising GC/MS, *n*-alkane concentrations may be overestimated (particularly for >C₁₉) and
41
42 352 thus may affect the calculated/reported ratios. According to previously published literature, one
43
44 353 possibility of compound X could be a dimethyl isomer belonging to the C_{*n*+1}-alkane.^{43,60} Dimethyl
45
46 354 isomers can be shifted by ~100 delta-Kovats (~1 carbon number),⁴³ whereas tri-methyl and tetra-
47
48 355 methyl have been reported to be shifted by ~150 and ~200 delta-Kovats.⁶⁰ The shift in delta-Kovats
49
50 356 for tetra-methyl alkanes identified in this study for example, 2,6,10,14-tetramethylhexadecane (C₂₀),
51
52 357 2,6,10,14-tetramethylheptadecane (C₂₁), 2,6,10,14-tetramethyloctadecane (C₂₂), 2,6,10,14-

1
2 358 tetramethylnonadecane (C₂₃, see Figure 4) etc., are in agreement with previous studies (~2 carbon
3
4 359 numbers)⁶⁰ and therefore are not likely candidates for compound X, due to the presence of the C_{n+1}-
5
6 360 alkane molecular ion in the mass spectra of the latter (see Figure 5).
7

8
9 361

10 362 *Monocyclic Alkanes and Aromatics*

11
12 363 *n*-Alkylcycloalkanes have been recognised as major constituents in sediments and crude oils.⁶¹
13
14 364 Worton et al.¹⁷ recently reported that primary organic aerosol (POA) emissions were dominated by
15
16 365 branched cycloalkanes with one or more rings and one or more branched alkyl side chain (>80% by
17
18 366 mass). The authors investigated the chemical composition of POA and its similarity to lubricating
19
20 367 oil. It was shown that cyclic and branched alkanes dominated the mass of POA (>60%), while
21
22 368 aromatics and *n*-alkanes accounted for <2 and <5%, respectively. The components of the base oil in
23
24 369 this study were predominantly cyclic and branched alkanes (based on TIC signal, see Supporting
25
26 370 Information), in agreement with Worton et al.¹⁷, but were not for the 5W30 motor oil, where
27
28 371 straight and branched chain alkanes were dominant.¹⁷ This is associated with the differences in
29
30 372 physical properties of the oils (*e.g.* viscosity and viscosity index) see Table S1. It has been
31
32 373 suggested that cycloalkanes are predominant compounds found in lubricating oils, owing to their
33
34 374 desirable properties making them efficient lubricants across a wide range of operating
35
36 375 temperatures.¹⁷
37
38
39
40
41

42 376

43
44 377 *n*-Alkylcyclopentanes and *n*-alkylcyclohexanes were identified by retention times of known
45
46 378 standard materials, mass spectral fragmentation at both ionisation energies and by comparison with
47
48 379 published data.¹⁷ *n*-Alkylcycloheptanes have received considerably less attention, probably due to
49
50 380 the difficulties associated in measuring these compounds with various sensitive detectors or owing
51
52 381 to them being present in trace amounts. Although *n*-alkylcycloheptanes were tentatively identified
53
54 382 they were represented by very small chromatographic peaks in comparison to *n*-alkylcyclohexanes
55
56 383 and *n*-alkylcyclopentanes. The selected mass fragment ions for *n*-alkylcyclopentanes, *n*-
57
58
59
60

1
2 384 alkylcyclohexanes and *n*-alkylcycloheptanes were m/z 68, m/z 82 and m/z 97 respectively. The
3
4 385 identified monocyclic alkanes ranged from C₁₂ to C₂₆ and their selected ion chromatograms are
5
6 386 shown in Figure 6. These could be readily identified by 70 eV ionisation.
7

8
9 387

10 388 Using 14eV fragmentation energy mass spectra, alkyl-methylcyclohexanes (m/z 82+97), alkyl-
11
12 389 dimethylcyclohexanes (m/z 82+111) and alkyl-trimethylcyclohexanes (m/z 82+125) were also
13
14
15 390 tentatively identified, (and confirmed by published retention indices for alkyl-methylcyclohexane)⁶²
16
17 391 using the summed selected ion chromatogram with their respective mass fragments as is illustrated
18
19 392 in Figure 6. The wealth of peaks identified signifies the complexity of the base motor oil with only
20
21 393 the monomethyl-alkylcycloalkanes being positively identified. Although it is not possible to
22
23 394 determine the positioning of the methyl group on the cyclohexane ring (*i.e.* ortho-, para- or meta-
24
25 395 position), it is possible to determine the positioning of branching on the alkyl chain as demonstrated
26
27 396 for the monomethyl alkanes. Three C₂₂H₄₄ isomers were differentiated by their mass spectral
28
29 397 differences (at 14eV), according to the positioning of the methyl group on the alkyl chain
30
31 398 (analogous to the monomethylalkanes) and the presence of a methyl group on the cyclohexane ring.
32
33 399 All three compounds, 1-methyl-4-pentadecylcyclohexane, (3-methylpentadecyl)cyclohexane and 1-
34
35 400 methyl-4-(3-methyltetradecyl)cyclohexane possessed their molecular ion of m/z 308, see Figure 7A
36
37 401 – 7C. 1-methyl-4-pentadecylcyclohexane and 1-methyl-4-(3-methyltetradecyl)cyclohexane showed
38
39 402 mass fragments of m/z 82+97, signifying a methyl-cyclohexane fragment. A relatively smaller peak
40
41 403 at m/z 97 was observed for (3-methylpentadecyl)cyclohexane signifying no methyl group on the
42
43 404 cyclohexane ring, and peaks m/z 111, 139, 169 and 197 were relatively large, corresponding to the
44
45 405 C-C bond breakage on either side of the methyl group on the alkyl chain, leading to mass fragments
46
47 406 C₈H₁₅, C₁₀H₁₉, C₁₂H₂₅ and C₁₄H₂₉, respectively, see Figure 7B. For compound 1-methyl-4-
48
49 407 pentadecylcyclohexane a relatively large peak at m/z 211 corresponding to the breaking for the alkyl
50
51 408 chain from the cyclohexane ring was observed, as well as peaks at m/z 111, 197, which most likely
52
53 409 signify a C₈H₁₅ ring fragment and a C₁₄H₂₉ aliphatic fragment, see Figure 7A. The compound 1-
54
55
56
57
58
59
60

1
2 410 methyl-4-(3-methyltetradecyl)cyclohexane, however, possessed relatively large peaks at m/z 125,
3
4 411 153, 155 and 183 which were significantly different to those of the other two isomers. These peaks
5
6 412 correspond to mass fragments C_9H_{17} , $C_{11}H_{21}$, $C_{11}H_{23}$ and $C_{13}H_{27}$, respectively, see Figure 7C.
7

8 413
9

10 414 Not all (methyl-alkyl)-cyclohexanes were positively identifiable in the motor oil, owing to the
11
12 415 intensities of the mass fragments m/z 82 and 97 being relatively large in comparison to the
13
14 416 individual bond cleavages along the alkyl chain. However, the methyl-alkylated cyclohexanes are
15
16 417 analogous in their positioning in the chromatograph to that of the n -alkanes and monomethylalkanes
17
18 418 (with higher polarity).¹⁷ Similarly to the alkyl-cycloalkanes, a series of alkyl-benzenes (m/z 92)
19
20 419 were positively identified ranging from tridecyl-benzene ($C_{19}H_{32}$) to eicosyl-benzene ($C_{26}H_{46}$).
21
22 420 However, surprisingly, no methyl-alkyl-benzenes were observed in either of the lubricating oils.
23
24 421 Alkyl-benzenes are often used as refrigeration lubricants and although they are present in mineral
25
26 422 oils, are found in small quantities. This is in agreement with the observed relative intensities of this
27
28 423 compound class in the chromatography, with respect to the cycloalkane derivatives which are more
29
30 424 abundant.¹⁷
31
32
33
34
35
36

37 425

38 426 *Bicyclic Alkanes*

39 427 Bicyclic aliphatics occupy the same section of the chromatogram as the monocyclic aliphatics,
40
41 428 albeit with higher polarity (larger in the y-axis of the chromatography). Alkyl-
42
43 429 decahydronaphthalene (or alkyl-decalin) isomers have received increasing interest as potential
44
45 430 petroleum biomarkers. Tran et al.⁶³ tentatively observed this compound class using GC×GC-ToF-
46
47 431 MS and stressed that the assignments were necessarily tentative owing to the lack of molecular
48
49 432 specificity in the interpretation of the mass spectra. Avila et al.⁶⁴ utilised GC×GC-ToF-MS and EI
50
51 433 coupled to Fourier transform cyclotron resonance mass spectrometry to analyse compounds found
52
53 434 in heavy gas oils, tentatively identifying alkyl-decalins. More recently, however, Wang et al.⁶⁵
54
55 435 characterised these compounds by synthesizing and investigating their mass spectra using GC×GC-
56
57
58
59
60

1
2 436 ToF-MS, reporting the mass spectra and retention times for decalins ranging from C₁ to C₁₁ in alkyl
3
4 437 chain length. Worton et al.¹⁷ reported that bicyclic alkanes with various branched alkyl chains
5
6 438 contributed about 30% of the observed mass in POA and lubricating oil. Up to 11 % of the observed
7
8 439 mass fraction of diesel fuel has also been attributed to bicyclic alkanes.⁴³ In this study, we identify
9
10 440 alkyl-decalins without the use of authentic standards owing to the characteristic mass spectra
11
12 441 obtained for these compounds, with the presence of their corresponding molecular ion; the mass
13
14 442 fragments of the respective observed isomers are in good agreement with the results published by
15
16 443 Wang et al.⁶⁵. Alkyl-decalins were identified by their characteristic mass fragment *m/z* 137 together
17
18 444 with the molecular ion (using 14eV spectra). The two closely eluting peaks correspond to the alkyl
19
20 445 chain at the 1- and 2- positions of the decahydronaphthalene structure but cannot be distinguished
21
22 446 from one another. No (methyl-alkyl)-decalins, however, were observed in either of the lubricating
23
24 447 oils, which is in contrast to other studies.¹⁷
25
26
27
28
29

448

449 *Tricyclic Aromatics and Alkanes*

32
33 450 Tricyclic aromatic compounds occupy the bottom half of the chromatographic section of the
34
35 451 contour plot illustrated in Figure 1B. In the base oil, these compounds were more abundant
36
37 452 (according to TIC signal, see Supporting Information) than the bicyclic aromatics/alicyclics but less
38
39 453 than the monocyclic aromatics/alicyclics. Selecting mass fragment *m/z* 191, two groups of tricyclic
40
41 454 aromatic compounds were identified; alkyl-anthracenes (or alkyl-phenanthrenes) and
42
43 455 alkyl-perhydroanthracenes (or alkyl-perhydrophenanthrenes). The mass spectra of alkyl-anthracenes
44
45 456 showed little fragmentation, possessing a dominant molecular ion (M⁺) and the *m/z* 191, signifying
46
47 457 the 3 fused benzene rings; whereas the observed mass fragments for alkyl-perhydroanthracenes
48
49 458 were *m/z* 55, 67, 81, 95, 109, 135 with the two most intensive peaks belonging to the M⁺ molecular
50
51 459 ion and *m/z* 191. These mass spectra were confirmed using the interpretation of Kiselev et al.⁶⁶, who
52
53 460 investigated the mass spectrometry of the stereoisomers of perhydroanthracene and
54
55 461 perhydrophenanthrene. Compounds from both tricyclic groups ranged from C₁ (methyl-) to C₁₂
56
57
58
59
60

1
2 462 (dodecyl-) in alkyl chain length, corresponding to molecular masses (M^+) of m/z 192 to m/z 360,
3
4 463 respectively. Very limited tricyclic aliphatics were found in either of the lubricating oils in contrast
5
6 464 with other studies.^{17,67} We attribute this to the differences in formulation of different oils, or to the
7
8 465 different methods used as highlighted in a recent study.⁶⁸

9
10 466

11 12 13 467 **CONCLUSION**

14
15 468 The primary reason why very little quantitative analysis has been conducted on hydrocarbons of
16
17 469 greater than C_{10} is that as a consequence of their huge diversity conventional GC is unable to
18
19 470 provide a separation and these compounds appear in the chromatogram as a large hump referred to
20
21 471 as UCM. The advent of two dimensional gas chromatography techniques in recent years has
22
23 472 provided a means to disaggregate the UCM hump providing the separation of thousands of
24
25 473 compounds which can be characterised on the basis of their mass spectra. This however has its
26
27 474 limitations as each organic compound, particularly in petroleum or environmental samples deriving
28
29 475 from petroleum, has a large number of isomers with rather similar mass spectra. The 70eV mass
30
31 476 spectra show minor differences in the fragmentation patterns which offer opportunities for
32
33 477 distinguishing isomeric species; however, a substantial limitation is that the mass spectra rarely
34
35 478 contain a signal from the molecular ion hence greatly complicating the identification process as
36
37 479 rather similar mass-spectra at 70eV can result from compounds of differing carbon chain length, as
38
39 480 demonstrated by Figure 4(A). Using a variable energy EI source operating at 14 eV has allowed the
40
41 481 generation of mass-spectra with clearly visible molecular ions. The combination of retention times
42
43 482 in two dimensions and mass-spectra at low and high ionisation energies confers unparalleled power
44
45 483 to identify specific isomers within the chromatograms.

46
47 484

48 485 **ACKNOWLEDGEMENTS**

49
50 486 The authors would like to thank L. McGregor, S. Smith and I. Biggs (Markes International) for
51
52 487 valuable discussions on aspects of GC×GC-ToF-MS. Financial support from the European
53
54
55
56
57
58
59
60

1
2 488 Research Council (ERC-2012-AdG, Proposal No. 320821) for the FASTER project is gratefully
3
4 489 acknowledged. The assistance of Zhirong Liang in aspects of data analysis and Soheil Zeraati-
5
6 490 Rezaei in the provision of oil samples is also acknowledged.
7

8 491
9

10 492 **Conflict of Interest Disclosure**

11
12 493 The authors declare no competing financial interest.
13
14

15 494
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

REFERENCES

- (1) Mader, P. P.; MacPhee, R. D.; Lofberg, R. T.; Larson, G. P. *Ind. Eng. Chem.* **1952**, *44*, 1352-1355.
- (2) Abelson, P. H. *Science* **1998**, *281*, 1609.
- (3) McDonald, T. A. *Chemosphere* **2002**, *46*, 745-755.
- (4) Kim, K. H.; Jahan, S. A.; Kabir, E.; Brown, R. J. C. *Environ. Int.* **2013**, *60*, 71-80.
- (5) Goldstein, A. H.; Galbally, I. E. *Environ. Sci. Technol.* **2007**, *41*, 1514-1521.
- (6) Ventura, G. T.; Kenig, F.; Reddy, C. M.; Frysinger, G. S.; Nelson, R. K.; Mooy, B. V.; Gaines, R. B. *Org. Geochem.* **2008**, *39*, 846-867.
- (7) Sutton, P. A.; Lewis, C. A.; Rowland, S. J. *Org. Geochem.* **2005**, *36*, 963-970.
- (8) Jacquot, F.; Doumenq, P.; Guiliano, M.; Munoz, D.; Guichard, J. R.; Mille, G. *Talanta* **1996**, *43*, 319-330.
- (9) Stout, S. A.; Wang, Z. *Oil Spill Environmental Forensics: Fingerprinting and Source Identification*, 1st Ed.; Academic Press: Burlington, 2007.
- (10) Schauer, J. J.; Kleeman, M. J.; Cass, G. R.; Simoneit, B. R. T. *Environ. Sci. Technol.* **1999**, *33*, 1578-1587.
- (11) Omar, N. J.; Abas, M. R.; Rahman, N.; Tahir, N.; Rushdi, A.; Simoneit, B. T. *Environ. Geol.* **2007**, *52*, 1485-1500.
- (12) Harrison, R. M.; Yin, J. J. *Environ. Monit.* **2010**, *12*, 1404-1414.
- (13) Alves, C.; Nunes, T.; Vicente, A.; Gonçalves, C.; Evtugina, M.; Marques, T.; Pio, C.; Bate-Epey, F. *Atmos. Res.* **2014**, *150*, 57-68.
- (14) Garcia de Oteyza, T.; Grimalt, J. O. *Environ. Pollut.* **2006**, *139*, 523-531.
- (15) Gough, M. A.; Rowland, S. J. *Nature* **1990**, *344*, 648-650.
- (16) Mao, D.; Weghe, H. V. D.; Lookman, R.; Vanermen, G.; Brucker, N. D.; Diels, L. *Fuel* **2009**, *88*, 312-318.
- (17) Worton, D. R.; Isaacman, G.; Gentner, D. R.; Dallmann, T. R.; Chan, A. W.; Ruehl, C.; Kirchstetter, T. W.; Wilson, K. R.; Harley, R. A.; Goldstein, A. H. *Environ. Sci. Technol.* **2014**, *48*, 3698-3706.
- (18) Lim, Y. B.; Ziemann, P. J. *Environ. Sci. Technol.* **2009**, *43*, 2328-2334.
- (19) Aumont, B.; Camredon, M.; Mouchel-Vallon, C.; La, S.; Ouzebidour, F.; Valorso, R.; Lee-Taylor, J.; Madronich, S. *Faraday Discuss.* **2013**, *165*, 105-122.
- (20) Gentner, D. R.; Isaacman, G.; Worton, D. R.; Chan, A. W. H.; Dallmann, T. R.; Davis, L.; Liu, S.; Day, D. A.; Russell, L. M.; Wilson, K. R.; Weber, R.; Guha, A.; Harley, R. A.; Goldstein, A. H. *Proc. Natl. Acad. Sci.* **2012**, *109*, 18318-18323.
- (21) Liu, Z.; Phillips, J. B. *J. Chromatogr. Sci.* **1991**, *29*, 227-231.
- (22) Phillips, J. B.; Xu, J. *J. Chromatogr. A* **1995**, *703*, 327-334.
- (23) Alam, M. S.; West, C. E.; Scarlett, A. G.; Rowland, S. J.; Harrison, R. M. *Atmos. Environ.* **2013**, *65*, 101-111.
- (24) Adahchour, M.; Beens, J.; Brinkman, U. A. T. *J. Chromatogr. A* **2008**, *1186*, 67-108.
- (25) Arsene, C.; Vione, D.; Grinberg, N.; Olariu, R. I. *J. Liq. Chromatogr. Relat. Technol.* **2011**, *34*, 1077-1111.
- (26) Hamilton, J. F. *J. Chromatogr. Sci.* **2010**, *48*, 274.
- (27) Zhang, W.; Zhu, S.; He, S.; Wang, Y. *J. Chromatogr. A* **2015**, *1380*, 162-170.
- (28) Hamilton, J. F.; Webb, P. J.; Lewis, A. C.; Hopkins, J. R.; Smith, S.; Davy, P. *Atmos. Chem. Phys.* **2004**, *4*, 1279-1290.
- (29) Maccoll, A. *Org. Mass Spectrom.* **1986**, *21*, 601-611.
- (30) Maccoll, A. *Org. Mass Spectrom.* **1982**, *17*, 1-9.
- (31) Munson, B. *Anal. Chem.* **1977**, *49*, 772A-775A.
- (32) Munson, M. S. B.; Field, F. H. *J. Am. Chem. Soc.* **1966**, *88*, 2621-2630.
- (33) Leinweber, P.; Schulten, H. R. *J. Anal. Appl. Pyrolysis* **1999**, *49*, 359-383.
- (34) Fenselau, C.; Wang, S. Y.; Brown, P. *Tetrahedron* **1970**, *26*, 5923-5927.
- (35) Hanley, L.; Zimmermann, R. *Anal. Chem.* **2009**, *81*, 4174-4182.
- (36) Boesl, U. *J. Phys. Chem.* **1991**, *95*, 2949-2962.
- (37) Hejazi, L.; Ebrahimi, D.; Hibbert, D. B.; Guilhaus, M. *Rapid Commun. Mass Spectrom.* **2009**, *23*, 2181-2189.
- (38) Hejazi, L.; Guilhaus, M.; Hibbert, D. B.; Ebrahimi, D. *Rapid Commun. Mass Spectrom.* **2015**, *29*, 91-99.
- (39) Wang, F. C.-Y.; Qian, K.; Green, L. A. *Anal. Chem.* **2005**, *77*, 2777-2785.
- (40) Mühlberger, F.; Hafner, K.; Kaesdorf, S.; Ferge, T.; Zimmermann, R. *Anal. Chem.* **2004**, *76*, 6753-6764.

- 1
2 548 (41) Bente, M.; Sklorz, M.; Streibel, T.; Zimmermann, R. *Anal. Chem.* **2008**, *80*, 8991-9004.
3 549 (42) Zimmermann, R.; Welthagen, W.; Gröger, T. *J. Chromatogr. A* **2008**, *1184*, 296-308.
4 550 (43) Isaacman, G.; Wilson, K. R.; Chan, A. W.; Worton, D. R.; Kimmel, J. R.; Nah, T.; Hohaus, T.; Gonin, M.;
5 551 Kroll, J. H.; Worsnop, D. R.; Goldstein, A. H. *Anal. Chem.* **2012**, *84*, 2335-2342.
6 552 (44) Chan, A. W. H.; Isaacman, G.; Wilson, K. R.; Worton, D. R.; Ruehl, C. R.; Nah, T.; Gentner, D. R.;
7 553 Dallmann, T. R.; Kirchstetter, T. W.; Harley, R. A.; Gilman, J. B.; Kuster, W. C.; de Gouw, J. A.; Offenberg, J.
8 554 H.; Kleindienst, T. E.; Lin, Y. H.; Rubitschun, C. L.; Surratt, J. D.; Hayes, P. L.; Jimenez, J. L.; Goldstein, A. H. *J.*
9 555 *Geophys. Res.: Atmos.* **2013**, *118*, 6783-6796.
10 556 (45) *Markes International Ltd - BenchTOF-Select Specification Sheet L-0008 (290615-4)* **2015**.
11 557 (46) *Markes International Ltd - Application Note 528*. **2014**.
12 558 (47) Krkošová, Ž.; Kubinec, R.; Addová, G.; Jurdáková, H.; Blaško, J.; Ostrovský, I.; Soják, L. *Pet. Coal* **2007**,
13 559 *49*, 51-62.
14 560 (48) Lu, H.; Peng, P. a.; Sun, Y. *Org. Geochem.* **2003**, *34*, 745-754.
15 561 (49) Huang, X.; Meyers, P. A.; Wu, W.; Jia, C.; Xie, S. *Org. Geochem.* **2011**, *42*, 156-165.
16 562 (50) Kenig, F.; Damsté, J. S. S.; Kock-van Dalen, A.; Rijpstra, W. I. C.; Huc, A. Y.; de Leeuw, J. W. *Geochim.*
17 563 *Cosmochim. Acta* **1995**, *59*, 2999-3015.
18 564 (51) Szabó, A. H.; Podolec, P.; Ferenczy, V.; Kubinec, R.; Blaško, J.; Soják, L.; Górová, R.; Addová, G.;
19 565 Ostrovský, I.; Višňovský, J. *J. Chromatogr. B* **2015**, *978*, 62-69.
20 566 (52) Hengerics Szabo, A.; Podolec, P.; Ferenczy, V.; Kubinec, R.; Blasko, J.; Sojak, L.; Gorova, R.; Addova, G.;
21 567 Ostrovsky, I.; Visnovsky, J.; Bierhanzl, V.; Cabala, R.; Amann, A. *J. Chromatogr. B: Anal. Technol. Biomed.*
22 568 *Life Sci.* **2015**, *978-979*, 62-69.
23 569 (53) Fenn, J.; Mann, M.; Meng, C.; Wong, S.; Whitehouse, C. *Science* **1989**, *246*, 64-71.
24 570 (54) Gilbert, A.; Yamada, K.; Yoshida, N. *Org. Geochem.* **2013**, *62*, 56-61.
25 571 (55) Kanzari, F.; Syakti, A. D.; Asia, L.; Malleret, L.; Mille, G.; Jamoussi, B.; Abderrabba, M.; Doumenq, P.
26 572 *Environ. Sci. Pollut. Res. Int.* **2012**, *19*, 559-576.
27 573 (56) Mille, G.; Asia, L.; Guiliano, M.; Malleret, L.; Doumenq, P. *Mar. Pollut. Bull.* **2007**, *54*, 566-575.
28 574 (57) Bray, E. E.; Evans, E. D. *Geochim Cosmochim. Acta* **1961**, *22*, 2-15.
29 575 (58) Rieley, G.; Collier, R. J.; Jones, D. M.; Eglinton, G.; Eakin, P. A.; Fallick, A. E. *Nature* **1991**, *352*, 425-427.
30 576 (59) Peters, K. E.; Walters, C. C.; Moldowan, J. M. *The biomarker guide: biomarkers and isotopes in the*
31 577 *environment and human history*, 2nd ed.; Cambridge University Press, 2005.
32 578 (60) Carlson, D. A.; Bernier, U. R.; Sutton, B. D. *J. Chem. Ecol.* **1998**, *24*, 1845-1865.
33 579 (61) Fowler, M. G.; Abolins, P.; Douglas, A. G. *Org. Geochem.* **1986**, *10*, 815-823.
34 580 (62) Wang, Z.; Fingas, M.; Li, K. *J. Chromatogr. Sci.* **1994**, *32*, 367-382.
35 581 (63) Tran, T. C.; Logan, G. A.; Grosjean, E.; Ryan, D.; Marriott, P. J. *Geochim. Cosmochim. Acta* **2010**, *74*,
36 582 6468-6484.
37 583 (64) Ávila, B. M. F.; Vaz, B. G.; Pereira, R.; Gomes, A. O.; Pereira, R. C. L.; Corilo, Y. E.; Simas, R. C.;
38 584 Nascimento, H. D. L.; Eberlin, M. N.; Azevedo, D. A. *Energy Fuels* **2012**, *26*, 5069-5079.
39 585 (65) Wang, H.; Zhang, S.; Weng, N.; Zhang, B.; Zhu, G.; Liu, L. *Analyst* **2015**, *140*, 4694-4701.
40 586 (66) Kiselev, A. V.; Kulikov, N. S.; Curthoys, G. *Chromatographia* **1984**, *18*, 297-304.
41 587 (67) Wingert, W. S. *Fuel* **1992**, *71*, 37-43.
42 588 (68) Worton, D. R.; Zhang, H.; Isaacman-VanWertz, G.; Chan, A. W.; Wilson, K. R.; Goldstein, A. H. *Environ.*
43 589 *Sci. Technol.* **2015**, *49*, 13130-13138.


590

591

592

593

594 FIGURES



595

596 **Figure 1. (A)** Total ion current chromatogram of a base lubricating oil examined by GC×GC-ToF-
597 MS illustrating improved chromatographic resolution by GC×GC compared with GC (white line in
598 background). **(B)** Contour plot of total ion current chromatogram of a base lubricating oil examined
599 by GC×GC-ToF-MS. Colours indicate the intensities of response (dark blue > blue > white).

600

601

602

603

604

605


606

607

608


609

610



611
612 **Figure 2.** Experimentally determined mass spectra of 2-methyleicosane (at 14eV, 70eV and NIST
613 reference), 3-methyleicosane (at 14eV, 70eV and NIST reference) and 4-, 5- and 6-methyleicosanes
614 (at 14eV and 70eV)

616



617

618 **Figure 3.** Contour plots of (A) 5W30 lubricating oil; and (B) mineral base oil; illustrating species
619 identified within two *n*-alkane species (C₂₀ – C₂₁).

620

621

622

623

624

625

626

627

628


629

630

631

632

633



634

635 **Figure 4.** (A) Experimentally determined mass spectra from base lubricating oil for (A) 2,6,10,14-
636 tetramethyl-nonadecane (upper panel) and 2,6,10,14,18-pentamethyl-nonadecane (lower panel) at
637 70 eV; (B) 2,6,10,14-tetramethyl-nonadecane (upper panel) and 2,6,10,14,18-pentamethyl-
638 nonadecane (lower panel) at 14 eV; and (C) reference spectra of 2,6,10,14-tetramethyl-nonadecane
639 (upper panel) and 2,6,10,14,18-pentamethyl-nonadecane (lower panel) from NIST library.

640

641

642

643


644

645


646

647

648



650
651 **Figure 5. (A)** Contour plot of synthetic motor oil 5W30 focusing on C₁₈ to C₂₂. Red circles indicate
652 *n*-alkanes, green circles indicate co-eluting species belonging to the isomer of the (*n*+1)-alkane.
653 Embedded mass spectra of C₂₀ *n*-alkane at 14eV and 70eV in top left and compound X in the
654 bottom right, **(B)** Contour plot of mineral base oil signifying no co-eluting interfering species.



662

663 **Figure 6.** Selected ion chromatogram for summed mass fragments m/z 68+82+97+111, illustrating
664 peaks identified for *n*-alkylcyclopentanes m/z 68 (black font); *n*-alkylcyclohexane m/z 82 (red font);
665 alkyl-methyl-cyclohexanes m/z 82+97 (yellow circles); alkyl-dimethyl-cyclohexanes m/z 82+111
666 (green circles); and alkyl-trimethyl-cyclohexanes m/z 82+125 (light blue circles).

667

668

669

670

671

672

673

674

675

676

677

678

679

680

681

682

683

684

685

686

687

688

689

690

691




Figure 7. 14eV ionisation energy mass spectra of (A) 1-methyl-4-pentadecylcyclohexane, (B) (3-methylpentadecyl)cyclohexane, and (C) 1-methyl-4-(3-methyltetradecyl)cyclohexane

692

693

694

695

696

697

698

699

700

701

702

703

704

705

706

707

708

709


710


711

712

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

713 For Table of Contents Only





1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60